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Please note

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Sleep and body measurements are for personal reference only. They do not replace a doctor's advice.

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
Designed at 20 cm × 20 cm. For a 6 cm booklet, print at 30%. Duplex, flip on the long edge. Trim inside the crop marks, fold the center line, and staple.

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This guide explains how to connect the Ailin Smart Ring to the Ailin Sleep app and sync sleep, heart rate, and blood oxygen. Follow the pages in order. It usually takes one or two minutes.

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Open https://xmianai.com in the phone browser, go to Download, and save the installer. Or open https://xmianai.com/xiaomian-latest.apk directly.

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Open the file and follow the prompts. If the phone asks for permission to install, tap Allow. Then open Ailin Sleep and sign in.

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Measurements are for your own reference. They are not for medical diagnosis or treatment.

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Make sure Bluetooth is on, permissions are allowed, and the ring is charged and beside the phone. Turn Bluetooth off and on, then tap Start scan again.

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Keep the ring nearby and powered. Turn Bluetooth off and on, then try again. Do not connect it from two phones at once.

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Connected, but no sleep from last night

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Wear it through the night, then sync in this app after you wake. Do not let another app read the ring first.

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Wear the ring on your finger, not on the charger. Relax your hand and wait about 10 to 60 seconds.

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If it still fails, call the number on the last page. Please say which phone you use and where you stopped.

设计页面 20 cm × 20 cm。若要 6 cm 成品，打印缩放选 30%。双面、长边翻转。裁切线以内为成品，中线对折后用书钉装订。
Designed at 20 cm × 20 cm. For a 6 cm booklet, print at 30%. Duplex, flip on the long edge. Trim inside the crop marks, fold the center line, and staple.

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Wear the ring and wait about 10 to 60 seconds. If the app asks you to take it off the charger, do that first.

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Do not measure while the ring is on the charger. Blood oxygen, temperature, and heart rate variability may be missing.

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Wear it through the night. In the morning, sync in Ailin Sleep before opening another app.

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If the battery stays on “--”, wait a few seconds or reconnect.

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Parts

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Transparent resin inner tray

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The smooth inner surface that rests on the finger.

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Two small windows on the inner band. Face them to the palm side.

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The round contact. The charging cable attaches here.

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The small light beside the charging port. Its color shows charging or a measurement.

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
There is no power button. When the magnetic cable attaches and the light flashes, the ring is on.

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
Red flashes while charging and turns off when full. Take the ring off the charger before pairing or measuring.

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Choose a comfortable size. The index, middle, or ring finger is recommended.

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Place the sensors, marked in orange, against the palm side of the finger.

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It should fit without sliding or feeling tight. Wear it before sleep and sync after you wake.

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The red light flashes while charging. Keep the cable attached until the ring is full.

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The light turns off when the ring is fully charged. You can take it off.

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The green light flashes during a heart rate test. Keep the ring on and wait about 10 to 60 seconds.

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The red light flashes during a blood oxygen test. Wear the ring on your finger, not on the charger.

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Open https://xmianai.com, go to Download, install the app, and sign in. Until you sign in, records cannot be saved to your account.

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Keep Bluetooth on until sync finishes. Do not use airplane mode. If the phone asks to find nearby devices or to use location, tap Allow. Location is used only to find the ring.

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Charge the ring, take it off the charger, and keep it within 1 meter of the phone.

设计页面 20 cm × 20 cm。若要 6 cm 成品，打印缩放选 30%。双面、长边翻转。裁切线以内为成品，中线对折后用书钉装订。
Designed at 20 cm × 20 cm. For a 6 cm booklet, print at 30%. Duplex, flip on the long edge. Trim inside the crop marks, fold the center line, and staple.

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Open Ailin Sleep and follow these labels in the app.

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Step 1. Open Devices. Tap Add a device, or the plus sign at the top right if you already have one.

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Step 2. Under Wearables, tap Ailin Smart Ring. Turn on Bluetooth, keep the ring close, and make sure it is powered on.

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Step 3. Tap Start scan and stay on this page. Scanning takes about ten seconds.

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Step 4. Under Devices found, tap Connect. Wait until it shows Connected and the battery level.

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If the ring does not appear, charge it, keep Bluetooth on, and scan again. Do not connect it to two phones at once.

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Tap Sync and keep the app open. You can then see the record on Devices, Home, or Reports.

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Sync in Ailin Sleep. Another app may clear records that have not been synced.

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After you disconnect, the ring stays in My Devices. Tap it to connect again.

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Keep one Bluetooth connection at a time. On a new phone, disconnect the old one, sign in, and pair again.

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